
VISICOIL™ MR

LINEAR FIDUCIAL MARKER

ENGLISH

Linear fiducial markers for soft tissue tumor localization

VISICOIL MR Pre-Loaded Fiducial Marker

Contains One (1) VISICOIL MR Fiducial Marker; One (1) Needle

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

DEVICE DESCRIPTION

VISICOIL MR is a coiled wire marker for soft tissue tumor location for enhanced accuracy in radiation therapy planning and delivery. The coil is made of high-purity platinum wire (99.95%) per ASTM B562, and is available in a variety of diameters (0.35 mm, 0.50 mm, 0.75 mm, and 1.00 mm) and lengths (0.50 cm, 1.00 cm, and 2.00 cm).

The device is pre-loaded in 304 Stainless Steel per ASTM A240, A666 needles, and are available in a variety of gauges (17g, 18g, 19g, 20g, 21g, 22g) and lengths (8 cm, 15 cm, 20 cm, and 30 cm). Device has a paraffin wax (Beeswax) (Bone wax) plug at the tip to keep the coil from coming out.

The multi-modality linear fiducial marker offers high visibility in ultrasound, computed tomography (CT) scan, X-ray, fluoroscopy, electric portal imaging devices (EPID) and Magnetic Resonance Imaging (MRI) modalities allowing for efficient MR/CT fusion and intensity-modulated radiation therapy (IMRT). The device is a passive device and is permanently implanted in the patient via percutaneous, endoscopic ultrasound (EUS), endobronchial ultrasound (EBUS), or bronchoscope.

Device is provided sterile for single use.

Pre-loaded VISICOIL MR is sterilized by ethylene oxide [STERILE]EO sterilization.

INDICATIONS FOR USE

The VISICOIL MR is indicated for use to radiographically mark soft tissue for future therapeutic procedures.

CONTRAINDICATIONS

None known.

POSSIBLE COMPLICATIONS

The VISICOIL MR may be placed in soft tissue, such as the tongue, tonsil, esophagus, lung, mediastinum, breast, liver, pancreas, prostate, bladder, cervix, vagina, and rectum.

The choice of a particular device must be carefully weighed against the patient's overall evaluation. Circumstances listed below may reduce the chance of a successful outcome.

- Allergy to the VISICOIL MR (Platinum) materials or foreign body sensitivity. Where material sensitivity is suspected, appropriate tests must be made prior to implantation.
- Migration - Fiducial markers may move away from where they were originally implanted. There must be sufficient time between the implantation of markers and treatment in order for the markers to stabilize.

Possible complications specific to the device may include:

- Rectal Implant: Coagulopathy or Anti-coagulants that cannot be stopped
- Trans-Rectal Implant for Prostate: Chronic Inflammatory Bowel Disease (Should go Trans-perineal)

Other complications associated with surgical procedure may include:

- Bleeding, Fever, Foreign body reactions necessitating medical intervention, Infection, Inflammation, Pain, Pancreatitis, Perforation, Peritonitis, Pneumothorax, Tumor Seeding.

WARNINGS

- The VISICOIL MR Fiducial Markers should slide freely within the introducer needle. Do not use excessive force to deploy.
- Pre-Loaded VISICOIL MR [STERILE]EO must not be re-sterilized.
- VISICOIL MR Fiducial Markers are for single use only. Never reuse any implant.
- Re-use of single-use devices creates a potential risk to the patient and/or the user. It may lead to contamination and/or impairment of functional capability. Contamination and/or the limited functionality of the device may lead to injury of the patient.
- Do not use the device after the expiration of the sterilization date or if the sterile packaging is damaged or has already been opened. Use the properly stored product by the date stated on the packaging.
- Prior to use read all instructions, precautions, and warnings. Failure to do so may result in severe patient injury.
- Do not alter the needle or any other kit/set component during insertion, use, or removal.

The stainless-steel needle is not intended for use with MRI guided implants.

CAUTIONS:

- This product may only be used by a qualified healthcare professional.
- Aseptic technique should be employed to avoid contamination of the needle portions that will be inserted through the skin and into tissue including the tips and shaft.
- Only introducer needles that have been tested to be MR Conditional or MR Safe should be used in the MR room for MR guided placement.
- The use of VISICOIL MR 1.00 mm is recommended only when using megavoltage / portal imaging systems. The visibility of the 0.75 mm, 0.50 mm, and 0.35 mm version are not optimal in this modality.

MR SAFETY INFORMATION

Non-clinical testing has demonstrated that the VISICOIL MR (Platinum) is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 Tesla and 3.0 Tesla, only.
- Maximum spatial field gradient of 3,000 gauss/cm or less.
- Maximum MR system reported the whole body averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes if scanning in the normal operating mode of operation for the MR system.

Under the scan conditions defined above, the VISICOIL™ MR Marker is expected to provide a maximum temperature rise of less than 1.6°C after 15 minutes of continuous scanning.

ARTIFACT TESTING

In non-clinical testing, the image artifact caused by the 1.00 mm 2 cm long VISICOIL MR extends approximately 8 mm from this implant when imaged with a gradient echo pulse sequence and a 3.0 T MR system.

DIRECTIONS FOR USE:

The information contained in this package insert is necessary but not sufficient for the use of this device. The device should be used by physicians trained on the procedure for which the device is intended. The techniques and procedures described do not represent all medically acceptable protocols, nor are they intended as a substitute for the physician's experience and judgment in treating any specific patient.

Pre-surgical Procedure

Position the patient appropriately. Scan patient with ultrasound, CT, X-ray, fluoroscopy, and/or MRI for anatomical landmarks and to ascertain needle puncture site, angle and expected depth to reach the targeted site(s).

NOTE: CONSULT INSTRUCTIONS FOR USE FROM ULTRASOUND, CT, X-RAY, FLUOROSCOPY, OR MRI IMAGING MACHINE MANUFACTURER.

- The procedure needs to be performed using aseptic technique to avoid contamination of the device. Disinfect the site thoroughly.

VISICOIL MR PRE-LOADED VERSIONS Preparation

- Prior to opening pouch, ensure coil has not been displaced into the pouch.
- Carefully remove preloaded needles from the sterile pouch using aseptic procedures. Ensure plastic protective tube is intact.
- Carefully remove the device assembly from the pouch. The marker should slide freely.

INSTRUCTIONS FOR USE

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- Inspect the implant needle and assembly for damage.
- If the plastic protective cover is dislodged, proceed with caution as the tip of the needle is sharp and delicate. Re-engage protector if necessary.
- Remove protective tube when ready for use.
- The patient is usually under general or local anesthesia.

Needle and Marker Location

- Communicate with Radiation Oncologist to identify optimal marker placement locations
- The needles are usually guided to their intended location under ultrasound, fluoroscopy or CT guidance. Align the tip of the needle with the desired marker location. Confirm needle/marker location.

Deploying the VISICOIL MR (Linear Deployment)

- Pull back on stylet handle exposing stylet lock and remove stylet lock.
- Carefully advance stylet into needle until there is slight resistance.
- The stylet is now up against the VISICOIL MR inside the needle.
- Securely hold stylet handle fixed in space (DO NOT ADVANCE).

- Securely grasp implantation needle HUB, carefully rotate and retract needle cannula while NOT moving the stylet handle. This will leave the VISICOIL MR marker(s) linear.
- A moment of resistance is expected as the coil breaks through the bone wax tip.

CAUTION: DO NOT INJECT THE VISICOIL MR if linear deployment is desired.

Advance the stylet if a collapsed deployment is preferred (Only with 0.35 mm size)

CAUTION: Ensure hub and stylet handle are completely in contact while gently rotating to ensure the marker separates from the delivery system.

Removing the Needle

- Once the VISICOIL MR marker is fully deployed, gently rotate needle assembly prior to removing.

CAUTION: Rotating the needle assembly ensures the residual bone wax disengages from the needle tip. Not rotating needle may result in undesired marker placement.

NOTE: If significant resistance is felt while the VISICOIL MR marker is inside of the needle, remove the entire system prior to deployment and discard.

NOTE: Discard needle in the appropriate sharps disposal container.

PACKAGING

The VISICOIL MR is a sterile device in the form of a coil. The device shipped sterile and packaged in a Mylar/Tyvek pouch as a single package. The device is sterilized with ethylene oxide sterilization. Handle with care.

HANDLING AND STORAGE

The VISICOIL MR should be stored under general office conditions. Protect the device from direct exposure to light, heat, and humidity.

Markers should be handled with caution as they are delicate devices and should be opened with care. Rough handling can cause damage to the device. Additional handling may result in damage. Uncontaminated defective markers and needles should be returned to RadioMed Corporation.

CUSTOMER COMPLAINTS
Communicate suspected deficiencies in product quality, identity, durability, reliability, safety, effectiveness and/or performance directly to RadioMed Corporation. Email: visicoil@iba-group.com Tel: For North America and the European Union, +1 (901) 432-7202.

When filing a complaint, please provide the component name(s), part number(s), lot number(s), your name and address, the nature of the complaint, and patient case number. Notify RadioMed Corporation immediately of an incident resulting in patient death or serious injury.

FURTHER INFORMATION
If further directions for the use of this system are needed, contact RadioMed Corporation Customer Service, email: visicoil@iba-group.com | Tel: +1 (901) 432-7202 or +1 (901) 805-5800.

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