VISICOIL™ MR

LINEAR FIDUCIAL MARKER

ENGLISH

Linear fiducial markers for soft tissue tumor localization

QTY: 1 Each

Contents: One (1) VISICOIL MR; One (1) Carrier Assembly

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician

DEVICE DESCRIPTION

VISICOIL MR is a coiled wire marker for soft tissue tumor location for enhanced accuracy in radiation therapy planning and delivery. The coil is made of high-purity platinum wire (99.95%) per ASTM B562, and is available in a variety of diameters (0.35 mm, 0.50 mm, 0.75 mm, and 1.10 mm) and lengths (0.50 cm, 1.00 cm, 2.00 cm, and 3.00cm).

The device is packaged on a 304 stainless steel carrier in a Mylar/Tyvek pouch.

The multi-modality linear fiducial marker offers high visibility in ultrasound, computed tomography (CT) scan, X-ray, fluoroscopy, electric portal imaging devices (EPID) and MRI image modalities allowing for efficient MR/CT fusion and intensity-modulated radiation therapy (IMRT).

The dense mass of the VISICOIL MR is recommended for improved imaging in all x-ray based imaging devices and other electromagnetic properties for enhanced visibility within Magnetic Resonance Imaging (MRI). The device is a passive device and is permanently implanted in the patient via percutaneous, endoscopic ultrasound (EUS), endobronchial ultrasound (EBUS), or bronchoscope.

The device is provided sterile for single use by gamma [STERILE] sterilization.

INDICATIONS FOR USE

The VISICOIL MR is indicated for use to radiographically mark soft tissue for future therapeutic procedures.

CONTRAINDICATIONS

None known.

POSSIBLE COMPLICATIONS

The VISICOIL MR may be placed in soft tissue, such as the tongue, tonsil, esophagus, lung, mediastinum, breast, liver, pancreas, prostate, bladder, cervix, vagina, and rectum.

The choice of a particular device must be carefully weighed against the patient's overall evaluation. Circumstances listed below may reduce the chance of a successful outcome.

- Allergy to the VISICOIL MR (Platinum) materials or foreign body sensitivity. Where material sensitivity is suspected, appropriate tests must be made prior to implantation.

- Migration - Fiducial markers may move away from where they were originally implanted. There must be sufficient time between the implantation of markers and treatment in order for the markers to stabilize.

Possible complications specific to the device may include:

- Rectal Implant: Coagulopathy or Anti-coagulants that cannot be stopped.
- Trans-Rectal Implant for Prostate: Chronic Inflammatory Bowel Disease (Should go Trans-perineal).

Other complications associated with surgical procedure may include:

Bleeding, Fever, Foreign body reactions necessitating medical intervention, Infection, Inflammation, Pain, Pancreatitis, Perforation, Peritonitis, Pneumothorax, Tumor Seeding.

WARNINGS

- VISICOIL MR [STERILE] must not be re-sterilized.

- Device for single use only. Never reuse any implant.

- Re-use of single-use devices creates a potential risk to the patient and/or the user. It may lead to contamination and/or impairment of functional capability. Contamination and/or the limited functionality of the device may lead to injury of the patient.

- Do not use the device after the expiration of the sterilization date or if the sterile packaging is damaged or has already been opened. Use the properly stored product by the date stated on the packaging.

- Prior to use read all instructions, precautions, and warnings. Failure to do so may result in severe patient injury.

- Do not alter the needle or any other kit/set component during insertion, use, or removal.

- Stainless-steel carrier has not been MRI Safety tested and is not intended for use with MRI guided implants.**

CAUTIONS:

- This product may only be used by a qualified healthcare professional.

- Aseptic technique should be employed to avoid contamination of the needle portions that will be inserted through the skin and into tissue including the tips and shaft.

- Only introducer needles that have been tested to be MR Conditional or MR Safe should be used in the MR room for MR guided placement.

- The use of VISICOIL MR 1.00 mm is only recommended when using megavoltage / portal imaging systems. The visibility of the 0.75 mm, 0.50 mm, and 0.35 mm version are not optimal in this modality.

MR SAFETY INFORMATION

Non-clinical testing has demonstrated that the VISICOIL MR Marker (Platinum) is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 Tesla and 3.0 Tesla, only.

- Maximum spatial field gradient of 3,000 gauss/cm or less.

- Maximum MR system reported the whole body averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes if scanning in the normal operating mode of operation for the MR system.

Under the scan conditions defined above, the VISICOIL MR Marker is expected to produce a maximum temperature rise of less than 1.6°C after 15 minutes of continuous scanning.

CAUTION: The stainless-steel carrier has not been MRI Safety tested and is not intended for use with MRI guided implants.

ARTIFACT TESTING

In non-clinical testing, the image artifact caused by the 1.00 mm 2 cm long VISICOIL MR Marker extends approximately 8 mm from this implant when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

DIRECTIONS FOR USE:

The information contained in this package insert is necessary but not sufficient for the use of this device. The device should be used by physicians trained on the procedure for which the device is intended. The techniques and procedures described do not represent all medically acceptable protocols, nor are they intended as a substitute for the physician's experience and judgment in treating any specific patient.

This instruction for use is focused on opening and loading device into a needle for implant into a patient. For detailed instruction on how to deploy, please refer to VISICOIL MR Pre-Loaded instructions for use.

CAUTION: Use needles with Flat or Ball Tipped stylet.

Step-By-Step Loading Guide Preparation

The Carrier Assembly is designed to aid loading the coil into a needle.

Hub Loading

- 1) Insert needle tip into sterile Bone Wax to block tip of needle.
- 2) Store needle in safe procedure and set aside within sterile field (This will be used later).

Figure 1: Removing the protective PE Tube

- 3) Remove and discard the Protective PE Tube. Hold Carrier vertical with Red Tip Protector upward. Firmly hold Carrier Hub and gently pull Protective PE Tube upward. (See Figure 1)

CAUTION: Keep tip of Carrier oriented upward as Marker is free to slide off. **CAUTION:** DO NOT contact the Marker with the inside of Protective Tube.

CAUTION: If tube contacts marker, inspect to ensure marker does not lift off with Tube.

- 4) Hold Needle horizontal and insert Carrier into Needle Hub until Carrier Hub rests on Needle Hub (See Figure 2).

CAUTION: If carrier must be inverted (Pointing downward) to insert into needle hub, use sterile, gloved index finger on coil to retain coil on carrier until inserted into needle hub.

Figure 2: Insert the carrier into the needle hub

- 5) Rotate Needle and Carrier vertical, so the Carrier is on top (See Figure 3).

- 6) Unscrew Stylet Wire Screw Cap and remove, allowing marker to fall into Needle (See Figure 3).

Figure 3: Removing the carrier delivery system

- 7) Remove and discard the entire Carrier Assembly and Screw Cap.

- 8) Re-insert and advance Stylet from Step 1, until slight resistance is present and re-apply Stylet Lock if available.

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until slight resistance is present and re-apply Stylet Lock if available.

CAUTION: Do not over-insert stylet as the Marker may be pushed out prematurely.

NOTE: To deploy device from a Needle, refer to Pre-Loaded Instructions for Use.

PACKAGING
The VISICOIL MR Markers are sterile devices in the form of a coil. The device is shipped sterile and packaged in a Mylar/Tyvek pouch as a single package. The device is sterilized with gamma sterilization. Handle with care.

HANDLING AND STORAGE
Device should be stored under general office conditions. Protect device from direct exposure to light, heat, and humidity.

Device should be handled with caution as they are delicate and should be opened with care.

Rough handling can cause damage to the device. Additional handling of the marker may result in damage. Uncontaminated defective product should be returned to RadioMed Corporation.

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