



Product Information

Multi-Modality Fiducial Markers for patient registration under CT and MR



Multi-Modality Fiducial Markers are self-adhesive skin markers which appear on CT and MR diagnostic imaging scans. The markers are used in image-guided surgery procedures as reference points during patient registration (alignment of 3D patient data with patient anatomy).

- for CT and MR
- soft material for patient comfort
- Center hole for patient marking

Item	Description	REF
Multi-Modality Fiducial Markers (60 pcs.)	• 60 Marker per box	MM3002-60
		

Our Multi-Modality Fiducial Markers stick out by:

All-purpose

The markers appear clearly on CT and MR diagnostic imaging scans

Adhesion

Best adhesion on skin - even with over-night patients

Functionality

2 mm center hole enables patient marking

Patient comfort

Soft hydrogel allows high patient comfort, even over-night