




IZI Medical

Osteo-Site® Murphy M1M™ - Side Bevel Special Size

The Osteo-Site® Murphy M1M Side Bevel is intended for vertebral body biopsy and infusion. The outer cannula of the needle is back-beveled to achieve improved purchase on the periosteum. Also, the side-bevel allows better direction of the needle to the desired location for more accurate results.



Key Features and Benefits:

- Provides longer lateral access to the bone¹
- Inner stylet and outer cannula are side beveled to match
- A blunt stainless steel obturator is included
- A 20 cm long radiolucent needle holder is included

References:

1. Single-Needle Lateral Sacroplasty Technique; Nicholson, P.J., et al., Spine, February 2019.

SET ORDER NUMBER	OUTER CANNULA GAUGE • LENGTH	GLOBAL PRODUCT NUMBER
DBBN-11-18.0-M1M	11G • 18cm	G27000



izimed.com

