



VISICOILTM

LINEAR FIDUCIAL MARKER




IZI Medical

Radiation Therapy Products

Products for Peak Radiation Therapy Performance.

Radiation oncologists, therapists, oncology nurses, and technologists like you understand the value of time and accuracy. That's why **IZI Medical Products®** ensures that our tools reduce procedural time and help you precisely target malignant tissue, creating a better overall experience for you, your team, and your patients.

We take pride in watching our concepts become products that transform and enhance your ability to perform the highest quality simulation and radiotherapy treatments. From simulation through treatment, our products ensure optimal patient stabilization and comfort while providing all the marking tools you need to better verify localization and reduce your number of fractions. Our innovation and dedication to accuracy are why so many healthcare professionals order our radiotherapy products by name.



Quality Care is our Primary Concern

IZI Medical Products is constantly striving for innovation by continually focusing on developing better products for you, and in the end for the benefit of your patients.

We understand the complex business of healthcare, and we know how important it is to balance quality, with efficiency and affordability.

As a leading innovator, manufacturer, and distributor of quality medical devices, we ensure that all of our products reduce procedural time and enhance precision. The result of which, is a better overall experience for you, your team and, most importantly, your patients. We take pride in watching our concepts take shape, transforming and enhancing the way clinical procedures are done. Our innovation and dedication to accuracy is why so many healthcare professionals order our products by name.

IZI Medical embraces the legacy of RadioMed Corporation's Visicoil fiducial marker portfolio, and looks forward to working together to improve patient care.

Visicoil™ Pre-Loaded Fiducial Markers

Linear Fiducial Markers for Soft Tissue Tumor Localization

The Pre-Loaded VISICOIL™ and Pre-Loaded VISICOIL^{MR} are available in various needle gauges with multiple marker sizes and lengths. This helps you choose the right implant method for your patient giving the most comfort and safety. The correct marker will optimally reduce common image distortions improving the accuracy of your treatment.

Product Benefits

- Wide range of needle sizes for physician preference
- Centimeter marks for depth placement
- Echogenic needle tip provides improved visualization
- Pre-Loaded, packaged sterile and ready for use
- Platinum = Improved MRI visibility

Product Specifications

- *Platinum Marker Diameter: 0.50mm, 0.75mm, 1.00mm
- Platinum Marker Length: 5mm, 10mm
- Gold Marker Diameter: 0.35mm, 0.50mm, 0.75mm, 1.10mm
- Gold Marker Length: 5mm, 10mm, 20mm, 30mm
- Needle Gauge: 17g, 18g, 19g, 20g
- Needle Length: 12cm, 20cm, 30cm

*Platinum Marker 0.35mm diameter available in the Interventional Radiology (SN) section



ORDER INFORMATION

Platinum

SET ORDER NUMBER	OUTER CANNULA GAUGE • LENGTH	VISICOIL™ DIAMETER • LENGTH
PRE-LOADED		
MR-050-005-PL	20G • 20cm	0.50mm • 5mm
MR-050-010-PL	20G • 20cm	0.50mm • 10mm
MR-075-005-PL	19G • 20cm	0.75mm • 5mm
MR-075-010-PL	19G • 20cm	0.75mm • 10mm
MR-100-005-PL	18G • 20cm	1.00mm • 5mm
MR-100-010-PL	18G • 20cm	1.00mm • 10mm

ORDER INFORMATION

Gold

SET ORDER NUMBER	OUTER CANNULA GAUGE • LENGTH	VISICOIL™ DIAMETER • LENGTH
PRE-LOADED		
VC-035-005-PL	19G • 20 cm	0.35mm • 5mm
VC-035-010-PL	19G • 20 cm	0.35mm • 10mm
VC-035-020-PL	19G • 20cm	0.35mm • 20mm
VC-035-030-PL	19G • 20cm	0.35mm • 30mm
VC-050-005-PL	18G • 20cm	0.50mm • 5mm
VC-050-010-PL	18G • 20cm	0.50mm • 10mm
VC-050-020-PL	18G • 20cm	0.50mm • 20mm

ORDER INFORMATION

Gold

SET ORDER NUMBER	OUTER CANNULA GAUGE • LENGTH	VISICOIL™ DIAMETER • LENGTH
PRE-LOADED		
VC-075-005-PL	18G • 20cm	0.75mm • 5mm
VC-075-010-PL	18G • 20cm	0.75mm • 10mm
VC-075-020-PL	18G • 20cm	0.75mm • 20mm
VC-075-030-PL	18G • 20cm	0.75mm • 30mm
VC-110-005-PL	17G • 20cm	1.10mm • 5mm
VC-110-010-PL	17G • 20cm	1.10mm • 10mm
VC-110-020-PL	17G • 20cm	1.10mm • 20mm
PRE-LOADED		
VC-050-005-PLSL	19G • 20cm	0.50mm • 5mm
VC-050-010-PLSL	19G • 20cm	0.50mm • 10mm
VC-050-020-PLSL	19G • 20cm	0.50mm • 20mm
VC-050-005-PL30	18G • 30cm	0.50mm • 5mm
VC-050-010-PL30	18G • 30cm	0.50mm • 10mm
VC-075-005-PL12	18G • 12 cm	0.75mm • 5mm
VC-075-010-PL12	18G • 12 cm	0.75mm • 10mm
VC-075-005-PL30	18G • 30cm	0.75mm • 5mm
VC-075-010-PL30	18G • 30cm	0.75mm • 10mm
VC-075-020-PL30	18G • 30cm	0.75mm • 20mm
VC-110-010-PL30	17G • 30cm	1.10mm • 10mm
VC-110-020-PL30	17G • 30cm	1.10mm • 20mm



Visicoil™ Interventional Radiology Fiducial Markers

Linear Fiducial Markers for Soft Tissue Tumor Localization

VISICOIL™ / VISICOIL^{MR} 0.50mm and 0.35mm markers are available pre-loaded in 21 and 22 gauge needles. Ideal for Interventional Radiology allowing for easier access with improved patient comfort and safety due to the smaller needle size.

Product Benefits

- Designed and approved for Interventional Radiology
- Reduced risk of pneumothorax
- 50% reduction in size vs. traditional gold fiducials
- Shorter needles ideal for standard bore CT
- Steerable needle aids in accurate marker placement
- Centimeter marks for depth placement

Product Specifications

- Platinum Marker Diameter: 0.35mm, 0.50mm
- Platinum Marker Length: 5mm, 10mm, 20mm
- Gold Marker Diameter: 0.35mm, 0.50mm
- Gold Marker Length: 5mm, 10mm, 20mm
- Needle Gauge: 21g, 22g
- Needle Length: 8cm, 15cm, 20cm



ORDER INFORMATION

Platinum

SET ORDER NUMBER	OUTER CANNULA GAUGE • LENGTH	VISICOIL™ DIAMETER • LENGTH
INTERVENTIONAL RADIOLOGY		
SNMR-035-005-2208	22G • 8cm	0.35mm • 5mm
SNMR-035-005-2215	22G • 15cm	0.35mm • 5mm
SNMR-035-005-2220	22G • 20cm	0.35mm • 5mm
SNMR-035-010-2208	22G • 8cm	0.35mm • 10mm
SNMR-035-010-2215	22G • 15cm	0.35mm • 10mm
SNMR-035-010-2220	22G • 20cm	0.35mm • 10mm
SNMR-035-020-2215	22G • 15cm	0.35mm • 20cm
SNMR-035-020-2220	22G • 20cm	0.35mm • 20mm
SNMR-050-005-2108	21G • 8cm	0.50mm • 5mm
SNMR-050-005-2115	21G • 15cm	0.50mm • 5mm
SNMR-050-005-2120	21G • 20cm	0.50mm • 5mm
SNMR-050-010-2108	21G • 8cm	0.50mm • 10mm
SNMR-050-010-2115	21G • 15cm	0.50mm • 10mm
SNMR-050-010-2120	21G • 20cm	0.50mm • 10mm



ORDER INFORMATION

Gold

SET ORDER NUMBER	OUTER CANNULA GAUGE • LENGTH	VISICOIL™ DIAMETER • LENGTH
INTERVENTIONAL RADIOLOGY		
SN-035-005-PL08	22G • 8cm	0.35mm • 5mm
SN-035-005-PL15	22G • 15cm	0.35mm • 5mm
SN-035-010-PL08	22G • 8cm	0.35mm • 10mm
SN-035-010-PL15	22G • 15cm	0.35mm • 10mm
SN-035-020-PL15	22G • 15cm	0.35mm • 20mm
SN-050-005-PL08	21G • 8cm	0.50mm • 5mm
SN-050-005-PL15	21G • 15cm	0.50mm • 5mm
SN-050-010-PL08	21G • 8cm	0.50mm • 10mm
SN-050-010-PL15	21G • 15cm	0.50mm • 10mm
SN-050-020-PL15	21G • 15cm	0.50mm • 20mm



Visicoil™ TwinLine Fiducial Markers*

Linear Fiducial Markers for Soft Tissue Tumor Localization

VISICOIL™ TwinLine is designed for image-guided robotic stereotactic radiosurgery. TwinLine stranded and non-stranded contain 2 VISICOIL™ markers with a specific length spacer, pre-loaded, sterile and ready to use. TwinLine consistent spacer gives reliable daily treatment set up while including the features of the traditional VISICOIL™. TwinLine requires less needles giving your patient more comfort and safety.

Product Benefits

- Designed for Cyberknife** image-guided robotic stereotactic radiosurgery
- Two markers, one needle
- Reliable consistent marker spacing
- Stranded / Non-stranded
- Faster, efficient daily set up and treatment
- Platinum offers improved MRI visibility

Product Specifications

- Marker Diameter: 0.50mm, 0.75mm
- Marker Length: 5mm
- Spacer Non-Stranded: 15mm, 20mm
- Spacer Stranded: 15mm, 20mm
- Needle Gauge: 18g, 19g, 20g, 21g
- Needle Length: 12cm, 20cm

*Visicoil TwinLine and Visicoil TwinLine MR devices are not available in all markets.

**Use of the Cyberknife* registered trademark by IZI Medical Products is not an endorsement by incorporation of the Visicoil TwinLine or its promotion and advertising. Cyberknife* is a trademark of Accuray Inc.



ORDER INFORMATION

Platinum

SET ORDER NUMBER	OUTER CANNULA GAUGE • LENGTH	VISICOIL™ DIAMETER • LENGTH	SPACER LENGTH	SLEEVE
TWINLINE MR *				
MRTL-050-41520	20G • 20cm	0.50mm • 5mm	15mm	
MRTL-050-41521	21G • 20cm	0.50mm • 5mm	15mm	
PTTL-050-11021	21G • 20cm	0.50mm • 5mm	10mm	
MRTLS-075-51518	18G • 20cm	0.75mm • 5mm	15mm	Yes
MRTLS-075-41519	19G • 20cm	0.75mm • 5mm	15mm	Yes

ORDER INFORMATION

Gold

SET ORDER NUMBER	OUTER CANNULA GAUGE • LENGTH	VISICOIL™ DIAMETER • LENGTH	SPACER LENGTH	SLEEVE
TWINLINE *				
TL-050-21520	20G • 20cm	0.50mm • 5mm	15mm	
TL-075-21518	18G • 20cm	0.75mm • 5mm	15mm	
TL-075-01518-12	18G • 12cm	0.75mm • 5mm	15mm	
TL-050-62020	20G • 20cm	0.50mm • 5mm	20mm	
TL-075-62018	18G • 20cm	0.75mm • 5mm	20mm	
TLS-075-31518	18G • 20cm	0.75mm • 5mm	15mm	Yes
TLS-075-31518-12	18G • 12cm	0.75mm • 5mm	15mm	Yes
TLS-075-72018	18G • 20cm	0.75mm • 5mm	20mm	Yes

*Visicoil TwinLine and Visicoil TwinLine MR devices are not available in all markets.



Visicoil™ Fiducial Markers

Linear Fiducial Markers for Soft Tissue Tumor Localization

IZI Medical Products offers VISICOIL™, the original flexible linear fiducial marker for use with today's most accurate image guided radiotherapy treatment plans. VISICOIL™/VISICOIL^{MR} are unique helically wound 99.95% pure gold or platinum fiducial markers. The coil design allows Visicoil™ to bond and flex with tissue providing reliable localization and set up during IGRT treatment. VISICOIL™ has been available since 2004 and has been used in hundreds of thousands of patients worldwide. Platinum having a higher magnetic susceptibility makes it ideal for use with MRI imaging.

Product Benefits

- Corrugated Surface = Enhanced ultrasound imaging
- Traditional Hollow Core = Minimal CT artifact
- Flexible Design = Excellent stability in target tissue
- Linear Design = Less needles and Faster implantation
- Platinum = Improved MRI visibility

Product Specifications

- Platinum Marker Diameter: 0.35mm, 0.50mm, 0.75mm, 1.00mm
- Platinum Marker Length: 5mm, 10mm
- Gold Marker Diameter: 0.35mm, 0.50mm, 0.75mm, 1.10mm
- Gold Marker Length: 5mm, 10mm, 20mm, 30mm



ORDER INFORMATION

Platinum

SET ORDER NUMBER	VISICOIL™ DIAMETER • LENGTH
LOOSE	
MRFP-035-005	0.35mm • 5mm
MRFP-035-010	0.35mm • 10mm
MRFP-050-005	0.50mm • 5mm
MRFP-050-010	0.50mm • 10mm
MRFP-075-005	0.75mm • 5mm
MRFP-075-010	0.75mm • 10mm
MRFP-100-005	1.00mm • 5mm
MRFP-100-010	1.00mm • 10mm

ORDER INFORMATION

Gold

SET ORDER NUMBER	VISICOIL™ DIAMETER • LENGTH
LOOSE	
FP-03505	0.35mm • 5mm
FP-0009	0.35mm • 10mm
FP-0010	0.35mm • 20mm
FP-0011	0.35mm • 30mm
FP-05005	0.50mm • 5mm
FP-05010	0.50mm • 10mm
FP-05020	0.50mm • 20mm
FP-0755	0.75mm • 5mm
FP-0015	0.75mm • 10mm
FP-0016	0.75mm • 20mm
FP-0017	0.75mm • 30mm
FP-11005	1.10mm • 5mm
FP-0057	1.10mm • 10mm
FP-0058	1.10mm • 20mm

About IZI Medical Products

Based in Owings Mills, Maryland, just outside Baltimore, IZI is a leading innovator, manufacturer, and distributor of quality medical devices that provides the continuity our healthcare partners need to succeed. Over the last three decades, IZI has built a diverse portfolio of products backed by strong intellectual property and currently sells to more than 1,000 domestic customers as well as internationally across 25 countries. IZI's portfolio of high-quality medical consumable devices are used in interventional radiology and oncology, radiation therapy, neuro-spine, and image guided surgery procedures.

Our mission is to improve each patient's quality of life by providing our healthcare professional partners a complete portfolio of clinically proven, high value, and innovative medical device solutions inspired by physicians, for physicians. With this mission in mind, our vision is to become a leader in diagnostic and therapeutic minimally invasive solutions to advance healthcare and improve patients' lives.

Please call for product samples, or visit izimed.com for more information.

- Radiation Therapy
- Diagnostic Radiology
- Image Guided Surgery
- Biopsy Localization
- Spine Solutions
- Breast Lesion Localization
- Biopsy
- Diagnostic Radiology



Primary Contact:

RadioMed Corporation
3150 Stage Post Drive / Suite 106
Bartlett, TN 38133 / USA

P / +1 901 432 7202
F / +1 901 432 7206
visicoil@izimed.com



IZI Medical

IZI Medical Products
5 Easter Court / Suite J
Owings Mills, MD 21117 / USA

P / +1 800 231 1499
P / +1 410 594 9403
IZIMed.com