

Breast Lesion Localization / X-Reidy


X-Reidy

Breast Lesion Localization Needle

Used for preoperative marking of nonpalpable breast lesions

- The X-shaped hook provides **excellent fixation** and is intended to **limit migration** before the procedure
- The 3 cm long stiff portion of the needle is **easily palpated** during the procedure
- The soft, flexible wire portion of the shaft is **comfortable** for the patient
- The reinforced wire is **easy to locate** during the procedure
- The reinforced wire is intended to **reduce the possibility** of accidentally cutting the wire


Supplied sterile in peel-open packages. Intended for one-time use only.



ORDER INFORMATION

X-Reidy - Standard Needles

SET ORDER NUMBER	OUTER CANNULA GAUGE • LENGTH	GLOBAL PROD. NO.
NEEDLE		
DXRBL-19.5-5.0-S	19.5G • 5 cm	G19617
DXRBL-19.5-9.0-S	19.5G • 9 cm	G19533
DXRBL-19.5-15.0-S	19.5G • 15 cm	G19589
SET		
DXRBL-18-5.0	18G • 5 cm	G19418
DXRBL-18-9.0	18G • 9 cm	G19380



Not all parts numbers shown may be approved for sale in all regulatory jurisdictions. Parts and configurations are subject to change without notification. Please consult with your Account Representative for details.