



IZI Medical

izimed.com

New VISICOIL^{MR}

Gold VISICOIL

NEW: VISICOIL MR 0.75MM X 0.6CM ON LEFT & 0.75MM X 0.6CM GOLD VISICOIL ON RIGHT

VISICOIL^{MR} PLATINUM
clearly visible in MRI

SIEMENS VERIO 3T 70CM OPEN BORE MRI, T2

VISICOILTMMR "new"

Linear fiducial marker for soft tissue tumor localization

Enhanced accuracy for advanced radiation therapy delivery

VISICOIL™ MR

Multi-modality linear fiducial marker for efficient MR/CT fusion and IGRT

SEE IT!


- Enhanced MR visibility = Faster MR registration and visualization
- Helical design offers excellent ultrasound imaging
- Visible under multiple imaging modalities to meet all your imaging needs
- Traditional Hollow Core = Minimal Artifact on both CT & MRI

TRUST IT!

- Stable marker localization from implant through treatment
- VISICOIL proven patented design
- Helical coil and flexible design fixates marker in target tissue
- Minimally invasive design allows for atraumatic placement

TREAT IT!

- Thin Needle Design = less pain, edema, bleeding, more comfort & safety
- Smaller needle path increases marker stability
- Linear Design = Less needles = Less measuring = Faster implantation
- Efficient, consistent, and reliable daily treatment set up



"Gold and Platinum have nearly identical high radiodensities, and are easily seen and tracked by X-ray based imaging. Gold is difficult to see on MRI, whereas Platinum has greater magnetic susceptibility, producing just the right signal void to be readily visible on MRI. Platinum is the ideal fiducial material."

DR. LEONARD AVRUCH

Department of Medical Imaging / University of Ottawa (Retired)

Available Sizes

0.35mm
0.50mm
0.75mm
1.00mm

0.5cm
1.0cm

Available preloaded in needles or loose

8cm needle
15cm needle
20cm needle



CONTACT US TODAY TO LEARN MORE

RadioMed Corporation

Bartlett, Tennessee 38133, USA

North America | T +1 901 432 7202

visicoil@izimed.com | izimed.com

